

BACKLIGHT KIT FOR 100 OZ BEAKERS



DESIGNED BY BAYER FOR EXCLUSIVE USE WITH CBP AND PSD TANKS

Backlight Kit-80013787

Designed to backlight 100 oz Top Mount Calibration Kit Beakers, mounted on CBP and PSD Tanks. Helps improve clarity when reading chemical level through heavy viscosity during calibration.

Features & Benefits

- Backlight mount (14.19" x 2.28")
- 12" led strip light
- On/off switch
- Standard 2-prong power cord (5ft)



Unlit

Backlit



GRADUATED BEAKER SIZE

- 100 Ounce - Height 22.28"

VALUE ADDED FEATURES

- Eye level mounting for easy viewing
- Easy field installation



Bayer SeedGrowth™

PRODUCT DETAILS: EASY INSTALLATION



1) Remove top wing nut

- Remove existing wing nut from the top of the Calibration Beaker



2) Loosen the lower nut

- Loosen the existing lower nut connected to the same stud as above



3) Insert the backlight bracket

- Slide the Backlight Kit bracket under the lower nut and over the top stud



4) Securely tighten in place

- Re-tighten the wing nut on the stud and the lower nut
- Connect the power cord to an external 115VAC power source



Bayer

Crop Science Division
1451 Dean Lakes Trail
Shakopee, Minnesota 55379 - USA
www.seedgrowth.bayer.com
Phone: +1 952 445 6868
Toll Free: 1 855 363 3152
E-mail: equipment@bayer.com



VALUE ADDED FEATURES

- Eye level mounting for easy viewing
- Easy field installation
- Light position accommodates any three existing beaker studs



Before using any crop protection product you should read and carefully follow directions, cautionary statements, and other information appearing on the product label.

Our technical information is based on extensive testing and is, to the best of our current knowledge, true and accurate but given without warranty as the conditions of use and storage are beyond our control. Descriptions and property data of the product do not contain any statement as to liability for possible damage.

SGRBACKLIGHTKIT01132016